

## PHOENIX TALKING SCALE MANUAL



Congratulations on purchasing a genuine My Weigh Phoenix Scale.

You now own a Personal Scale of proven technical design and highest quality.

We hope that you will enjoy the use of your new Personal Digital Scale.

You can visit us at [www.Myweigh.com](http://www.Myweigh.com)

### Function

Measures body weight and with voice reading.

### Technical specifications:

Weighing capacity:	150kg / 330lb
Minimum weighing limit:	2kg / 4.4 lb
Weighing graduation:	0.1kg / 0.2 lb
Weighing deviation:	$\pm (1\%+1d)$
Optimum temperature range:	20°C $\pm$ 15°C

### Operation instructions:

1. Inserting batteries:  
Insert a 9V battery into the battery compartment on the bottom (underside) of the scale.
2. Taking measurement
  - a. Step onto the measurement platform. The scale will say "Hello". Now step off the platform and wait a few moments
  - b. Once the platform is stable, it'll say "I'm ready" and show "0.0 lb" (or 0.0 kg). Then you may step onto the measurement platform again. Your weight will be displayed and read out once the reading is steady.
  - c. You may step off of the platform and the reading will hold for 10 secs. After that it'll say "Good-bye" and turn off automatically.
3. Power off  
The power will shut down automatically after no operation is detected for about 10 seconds.
4. Switching on/off the talking function  
Use the switch on the bottom of the platform to switch on/off the speaker function.
5. Unit conversion  
Turn the Scale ON, "0.0" will show on the display" press the kg/lb button on the front of the scale (under the display) to select the desired unit (LB or KG)
6. Important notes
  - \*\*\*If the Scale ever malfunctions please try installing a new battery, sometimes a low battery or defective battery may cause inaccurate readings or other errors
  - When the battery runs very low, the warning icon will appear on the display. When this happens, be sure to replace the battery immediately on the underside of the scale.
  - When over-loaded, "FFFF" will be displayed.
  - Use the unit on hard, flat, dry and secure surface only.
  - Place the unit in an area free from high temperature or high humidity. To clean, use dampened cloth and mild detergent.
  - Avoid excessive impact or vibration to the unit.
  - If you do not intend to use this unit for a long period of time, remove the battery before storage.
  - Do not attempt to disassemble the unit, otherwise the warranty will become void.
  - The unit is intended for home and personal use only (not legal for commerce).
  - This scale must only be used on a hard surface such as tile, wood linoleum or short carpet. This and most bath scales will not perform properly on thick-carpeted surfaces.



141-PHX-IM-03-01